



Laboratory diagnosis of fungal infections

Professor Ruoyu Li

Professor, Department of Dermatology

Peking University First Hospital

Peking University Research Center for Medical Mycology

Beijing, China



MMTN Conference 2019, April 18~19, Nanchang, China

Laboratory Diagnosis of Fungal Infection

Ruoyu Li, MD.

**Research Center for Medical Mycology, Peking University,
Beijing, China**

First Hospital, Research Center for Medical Mycology, Peking University



Introduction

First Hospital, Research Center for Medical Mycology, Peking University



Clinically important mycoses

- Skin and mucosal infections: ~billion
- Chronic severe (including NTDs) :
 - Chronic pulmonary aspergillosis
 - Mycetoma (NTD)
 - Chromoblastomycosis (NTD)
 - Sporotrichosis (NTD)
- Acute invasive
 - Candidiasis
 - Aspergillosis
 - Cryptococcosis
 - Mucomycosis
 - Pneumocystis jirovecii pneumonia
 - Teleromycosis marneffeii

First Hospital, Research Center for Medical Mycology, Peking University

Prevalence of fungal diseases-estimate precision

>150 million people have serious fungal diseases

- ~3,000,000 cases of chronic pulmonary aspergillosis
- ~223,100 cases of cryptococcal meningitis complicating HIV/AIDS
- ~700,000 cases of invasive candidiasis
- ~500,000 cases of PJP
- ~250,000 cases of invasive aspergillosis

Number of deaths from fungal disease:

>1.6million~tuberculosis

Bongomin F, Gago S, Oladele RO, et al. Global and multi-national prevalence of fungal diseases—estimate precision. J Fungi, 2017, 3(4): E57

First Hospital, Research Center for Medical Mycology, Peking University





Main determinants on incidence and prevalence of fungal disease

- **Socio-economic, geo-ecological characteristics ;**
- **Increasing number of at-risk populations :**
 - AIDS
 - Tuberculosis
 - COPD
 - Asthma
 - Cancers
 - Organ transplantation
 - Corticosteroid therapy

Bongomin F, Gago S, Oladele RO, et al. Global and multi-national prevalence of fungal diseases—estimate precision. J Fungi, 2017, 3(4): E57

First Hospital, Research Center for Medical Mycology, Peking University





Microbiological evidence-confirmed diagnosis

Microscopic exam.

- ☐ Sterile area ☐ Pathology

Culture:

- Positive in sterile area :
- ☐ CSF ☐ Tissues ☐ Extract
- ☐ Blood (Yeast, *Fusarium* etc., but not *Aspergillus*)

Antigen detection: Cryptococcus antigen

First Hospital, Research Center for Medical Mycology, Peking University



Microbiological evidence-Probable Diagnosis

Microscopic:

- Filamentous fungi positive in either:
- ☐ Sputum ☐ BAL ☐ Broncho-brush ☐ Sinus extract


Culture:

- Filamentous fungi positive in either :
- ☐ Sputum ☐ BAL ☐ Broncho-brush ☐ Sinus extract
- Antigen detection:
- G test (Serum)
- GM test (Plasma, Serum, BALF, CSF)

First Hospital, Research Center for Medical Mycology, Peking University



Diagnosis of Fungal Infection



First Hospital, Research Center for Medical Mycology, Peking University




Mycological Diagnosis

First Hospital, Research Center for Medical Mycology, Peking University



Direct Microscopic Examination

- KOH
- NS
- Calcofluor white (Fluorescent)
- India (Chinese) ink
- Gram
- Giemsa
- Histopathology



First Hospital, Research Center for Medical Mycology, Peking University





Materials

- | | |
|---------------------------------|-----------------|
| • Sputum | • Urine |
| • Bronchoalveolar lavage | • Fluids : |
| • Blood cultures | – Peritoneal |
| • Bone marrow | – Pleural |
| • Tissues and surgical material | – Pericardic |
| • Cerebrospinal fluid | – Ascitic fluid |
| | – Joint fluids |

First Hospital, Research Center for Medical Mycology, Peking University




Direct Examination-Yeast cell




First Hospital, Research Center for Medical Mycology, Peking University




Teleromycosis *marneffe*



intracellular yeast cell




First Hospital, Research Center for Medical Mycology, Peking University






Histopathology





Candidiasis



First Hospital, Research Center




Cryptococcosis-direct examination



First Hospital, Research Center for Medical Mycology, Peking University




Cryptococcosis-histopathology





First Hospital, Research Center for Medical Mycology, Peking University



Aspergillosis-direct examination





First Hospital, Research Center for Medical Mycology, Peking University





Mucomycosis-histopathology



First Hospital, Research Center for Medical Mycology, Peking University



Dermatophytosis-direct examination



First Hospital, Research Center for Medical Mycology, Peking University




CFW (Florescent) Staining

- Binding with polysaccharides on the chitin ring to develop fluorescence
- Combined with KOH dissolution
- Easier to find fungal elements
- Show the fungal structure clearly

First Hospital, Research Center for Medical Mycology, Peking University




Calcofluor white stain of *Malassezia*



First Hospital, Research Center for Medical Mycology, Peking University




Calcofluor white stain of dermatophytes



First Hospital, Research Center for Medical Mycology, Peking University




CFW staining of *Candida* species



First Hospital, Research Center for Medical Mycology, Peking University




CFW staining of *Aspergillus* species



First Hospital, Research Center for Medical Mycology, Peking University



CFW staining of *Mucor* species



First Hospital, Research Center for Medical Mycology, Peking University



Summary of direct examination

- Yeast cells: yeast;
- Yeast cells and pseudohyphae: *Candida* spp.;
- Yeast cells with capsules: *Cryptococcus* spp.;
- Transparent septate hyphae, about 45° branches: *Aspergillus* spp.;
- Transparent, non-septate hyphae, wide, about 90° branches: *Mucor* spp.;
- Brown or black hyphae or conidia: dematiaceous fungi

First Hospital, Research Center for Medical Mycology, Peking University



Different culture media



First Hospital, Rese

logy, Peking University

Colony morphology under different media



First Hospital, Research Center for Medical Mycology, Peking University



Explanation of the Culture Results

- *Cryptococcus* spp., *P.marneffe*: confirmed diagnosis;
- *Candida* spp., *Aspergillus* spp.: combine with clinical and other reference;
- Sterile area, blood and CSF: confirmed diagnosis;
- Pus, sputum or urine: carefully explanation, repeat the culture if necessary;
- Combination the direct exam and culture result is very important.

First Hospital, Research Center for Medical Mycology, Peking University




Which isolates need to identify into species level


- Isolates from the sterile area, blood, CSF, body fluid, BAL etc.;
- Isolates from ICU, burn pts., organ transplantation pts.,
- Exact ID is necessary to help the selection

First Hospital, Research Center for Medical Mycology, Peking University



Fungal ID-Differentiate yeast or mold





Principle of Yeast identification

- "Yeast" is not a formal taxa, but is a widespread form of growth in ascomycetes and basidiomycetes;
- Identification needs to combine morphological, physiological and biochemical characteristics as well as molecular and mass spectrometry characteristics;

First Hospital, Research Center for Medical Mycology, Peking University



Identification of molds

- Most molds can be identified by culture;
 - colony form, surface color, and growth rate
 - culture medium, temperature, etc.
- Transmitted culture to low nutrient medium to promote sporulation
 - PDA and CMA
- The characteristic of conidiogenesis under the microscope was observed by the technique of slide culture



First Hospital, Research Center for Medical Mycology, Peking University




Culture ID-Morphology

- The ID were determined by the morphology of colonies, microscopic characteristic of fungal sporulation;
- Try to identify under the microscope;
- Need to see sporulation (needle picking, tape or slide culture)
- Stimulate sporulation media: PDA and CMA

First Hospital, Research Center for Medical Mycology, Peking University






Culture ID-biochemical method

- CHROMagar Candida(48hr)
- API 20C AUX(48-72hr)
- RapID Yeast Plus System(4-5hr)
- Automatic system
 - ID 32C strip system(24-48hr)
 - Vitek Yeast Biochemical Card system
 - Vitek 2 ID-YST card system (24hr)
 - Quantum II
 - Biolog YT MicroPlate system (Biolog, USA)


First Hospital, Research Center for Medical Mycology, Peking University



Biochemical ID



API 20C AUX




ID 32C

First Hospital, Research Center for Medical Mycology, Peking University



CHROMagar Candida



First Hospital, Research Center for Medical Mycology, Peking University



Vitek 2 ID-YST card system



First Hospital, Research Center for Medical Mycology, Peking University



Culture ID—Molecular Way

- PCR-EIA
- RLB, AFLP, SSCP, RAPD
- Gene Chips
- PNA-FISH (fluorescein-labeled peptide nucleic acid fluorescent in-situ hybridization)
 - 26S rRNA, *Candida albicans*, *C.glabrata*
 - sensitivity :100%, specificity: 100%

First Hospital, Research Center for Medical Mycology, Peking University



Common gene target for molecular ID

- rRNA(rDNA)gene
 - Small Subunit (SSU) rDNA (18S)
 - Large Subunit (LSU) rDNA (26S, 28S)
 - Internal Transcribed Spacer (ITS)
 - Intergenic Spacer (IGS)





Others:

- elongation factor
- β -tubulin
- Calmodulin

First Hospital, Research Center for Medical Mycology, Peking University

Classification of pathogenic fungi based on 18S-rRNA





MALDI-TOF MS and other techniques


	Methods	Time	Cost (Yuan/sample)	Difficulties	Throughout
Conventional method	Phenotypic	7-14d	200	Complicated	low
	rDNA	1-2d	80	Complicated	low
	DNA-DNA hybridization	1-2d	1000	Complicated	5/day
	Lipid acid analysis	1-2h	1000	Special training	10/day
	MALDI-TOF MS	3-5min	4	Simple, automatically	400/day

“Proteomic phenotyping is revolutionizing diagnostic mycology as fully reflecting species/morph varieties but often overcoming taxonomic hindrance.”

Chierico FD *et al. J of Proteomics* 2012;75:3314-30.



MALDI-TOF ID of *Candida* spp.



First Hospital, Research Center for Medical Mycology, Peking University




Non-culture diagnosis

First Hospital, Research Center for Medical Mycology, Peking University



Non-culture Diagnostic Methods



First Hospital, Research Center for Medical Mycology, Peking University

Comparison of Different Methods

	G test	GM test	LA test
Pathogen	Can. Asp. Pneumocystis, etc.	Asp.	Crypto.
Methods	Agglutination	ELISA	Latax
Samples	plasma	serum	serum、CSF
Time	4hrs.	6hrs.	30mins.
Sensitivity	60-90%	71%	90-100%
Specificity	60-90%	89%	≈100%

First Hospital, Research Center for Medical Mycology, Peking University



The advantage of serology diagnosis

- Fast report
- Increased sensitivity
- Easy to get sample
- The examined makers are usually diagnostic
- Could be used for the evaluation of severity and treatment effect


First Hospital, Research Center for Medical Mycology, Peking University



Cryptococcus antigen detection -agglutination test



First Hospital, Research Center for Medical Mycology, Peking University



β -1,3-D-Glucan: Indication

- **Not only for *Candida* infection**
 - *Aspergillus*、*Fusarium* and other molds (Miyazaki, J Clin Micro 33:3115, '95)
 - *Cryptococcus* negative (a-glucan in CW)
- **Positive result could not confirm the fungal species**
- **May be reduced by glucanase within 5 days, the monitor period should be 2/week**

First Hospital, Research Center for Medical Mycology, Peking University



GM test: EIA (Platelia)


Detecte galactomannan, (GM)

3. Detection – avidin-enzyme + substrate

2. Detecting Ab – anti-galactomannan
[biotin conjugated]

Ag – blood, urine, BAL

1. Capture antibody –
monoclonal anti - galactomannan



Specificity = >90% sensitivity = 75%




Summary of GM test

- Screening for GM = aid for the early (pre-emptive?) diagnosis of IA in prolonged neutropenic patients:
 - Cutoff = 0.5
 - Frequency: 2x/week
 - Caveats
- Serum、BALF、CSF samples
- GM + CT: high predictive value
- GM + β -D-glucan/PCR: increased specificity and sensitivity

First Hospital, Research Center for Medical Mycology, Peking University



Real-time PCR for detecting *Candida* spp.



Design

- Species selection
- DNA target selection
- Primers design
- Species/groups specific probe design

Set up of multiplex real-time PCR

Clinical material

- DNA extraction method

First Hospital, Research Center for Medical Mycology, Peking University



Thank you for your attention !

